

Product Data Sheet

Product Name: Galiximab

Cat. No.: GC74483

Chemical Properties

Cas. No. 357613-77-5

Formula

M.Wt

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Galiximab

Background

Galiximab (IDEC 114) is a primatized immunoglobulin G1 (IgG1) lambda monoclonal antibody directed against the CD80 antigen and blocks CD80–CD28 binding. Galiximab has variable regions are primatized (cynomologous monkeys), and the constant regions are human. Galiximab can be used in research of B-cell lymphoma.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Smith SM, et, al. The anti-CD80 primatized monoclonal antibody, galiximab, is well-tolerated but has limited activity in relapsed Hodgkin lymphoma: Cancer and Leukemia Group B 50602 (Alliance). Leuk Lymphoma. 2013 Jul;54(7):1405-10.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA