

Product Data Sheet

Product Name: SM1-71-R

Cat. No.: GC80867

Chemical Properties

Cas. No. 2088180-03-2

Formula $C_{24}H_{28}ClN_7O$

M.Wt 465.98


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

SM1-71 -R is a reversible analog of SM1-71. SM1-71 -R lacks the acrylamide warhead and therefore cannot form a covalent bond with the kinase. SM1-71 -R can be used as a control compound for SM1-71 [1].

In Vitro, SM1-71-R (1 µM; 2 h) results in inhibition of phosphorylation of p-AKTS473 in H23-KRASG12C cells [1].

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Rao S, et al. A multitargeted probe-based strategy to identify signaling vulnerabilities in cancers. J Biol Chem. 2019 May 24;294(21):8664-8673.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA