

Product Data Sheet

Product Name: TY 52156

Cat. No.: GC13720

Chemical Properties

Cas. No. 934369-14-9

Chemical Name N-(4-chlorophenyl)-3,3-dimethyl-2-oxobutanimidic acid, 2-(4-chlorophenyl)hydrazide

SMILES CC(C)(C)C(/C(NNC1=CC=C(Cl)C=C1)=N/C2=CC=C(Cl)C=C2)=O


Formula $C_{18}H_{19}Cl_2N_3O$ M.Wt 364.3

Solubility $\leq 15\text{mg/ml}$ in ethanol; 25mg/ml in DMSO; 30mg/ml in dimethyl formamide Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

TY 52156 (1-(4-chlorophenylhydrazono)-1-(4-chlorophenylamino)-3,3-dimethyl-2-butanone) is a novel S1P3 receptor antagonist [1].

Sphingosine 1-phosphate (S1P) is a bioactive lysophospholipid mediator mainly released from

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

activated platelets and implicated in many biological responses, such as angiogenesis, vascular development, and cardiovascular function [1].

In vitro: TY-52156 preferentially inhibited the S1P-induced increase in $[Ca^{2+}]_i$ in S1P3-CHO cells. TY-52156 competitively inhibited the dose-dependent $[Ca^{2+}]_i$ increase elicited by S1P in S1P3-CHO with the K_i value of ~ 110 nM for S1P3 receptor. TY-52156 showed submicromolar potency and a high degree of selectivity for S1P3 receptor. TY-52156 (10 μ M) inhibited 24 GPCRs and three ion channels by $< 30\%$. S1P dose-dependently decreased CF (Coronary Flow) in perfused rat heart. S1P dose-dependently induced vasoconstriction in isolated canine cerebral arteries. In HCASMCs, TY-52156 inhibited S1P-induced Rho activation.

In vivo: in SD rats, oral bioavailability of TY-52156 was $\sim 70.9\%$. In rats, pretreatment with TY-52156 significantly attenuated FTY-720-induced bradycardia, a broad agonist of S1P receptors. Pretreatment with TY-52156 prevented the FTY-720-P-induced increase in $[Ca^{2+}]_i$ in a dose-dependent manner.

Reference:

[1] Murakami A, Takasugi H, Ohnuma S, et al. Sphingosine 1-phosphate (S1P) regulates vascular contraction via S1P3 receptor: investigation based on a new S1P3 receptor antagonist[J]. Molecular pharmacology, 2010, 77(4): 704-713.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA