

Product Data Sheet

Product Name: E3 ligase Ligand 79

Cat. No.: GC80959

Chemical Properties

Cas. No. 2152672-96-1

Formula $C_{13}H_{14}N_4O_2$


M.Wt 258.28

Solubility Storage Store at 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

E3 ligase Ligand 79 is a CRBN E3 ligase ligand that can be used for synthesis of PROTACs [1].

References:

[1]. Zhiyong Yu. Bifunctional compounds containing 2,5-substituted pyrimidine derivatives for degrading cyclin-dependent kinase 2 and/or cyclin-dependent kinase 4 via ubiquitin proteasome pathway. WO2025240536A1. 2025-11-20.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA