

Product Data Sheet

Product Name: Filapixant

Cat. No.: GC63444

Chemical Properties

Cas. No. 1948232-63-0


Formula $C_{24}H_{26}F_3N_5O_3S$ M.Wt 521.56

Solubility DMSO : 146.67 mg/mL (281.21 mM; Storage Store at -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Filapixant is a purinoreceptor antagonist extracted from patent WO2016091776A1, example 348. Filapixant is the active reference substance of Eliapixant[1].

[1]. Davenport AJ, et, al. 1,3-thiazol-2-yl substituted benzamides. WO2016091776A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA