

Product Data Sheet

Product Name: Fmoc-Cys(SASRIN™ resin)-OH (200-400 mesh, 0.3-0.6 mmol/g)

Cat. No.: GA21562

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

This resin has been developed to allow C-terminal modification of cysteine-containing peptides. Fmoc-protected cysteine is linked to SASRIN™-resin by its thiol moiety, which thus is already protected. After derivatization of the carboxylic acid (

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA