

## Product Data Sheet

---

Product Name: Aluminum sulfate hydrate

Cat. No.: GB56948

### Chemical Properties

Cas. No. 17927-65-0

Formula  $\text{Al}_2(\text{SO}_4)_3 \cdot x\text{H}_2\text{O}$

M.Wt 342.15 (anhydrous basis)


Solubility

Storage Store at room temperature, keep dry and cool

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA