

Product Data Sheet

Product Name: H-Phenylalaninol

Cat. No.: GA10497

Chemical Properties

Cas. No. 3182-95-4

Chemical Name (2S)-2-amino-3-phenylpropan-1-ol

SMILES OC[C@H](N)CC1=CC=CC=C1

Formula $C_9H_{13}NO$

M.Wt 151.2


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

H-Phenylalaninol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA