

Product Data Sheet

Product Name: Bafisontamab

Cat. No.: GC74510

Chemical Properties

Cas. No. 2437210-79-0

Formula M.Wt 237.63 kDa

Solubility Storage Store at -80°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Bafisontamab

Background

Bafisontamab (EMB-01) is a bispecific antibody targeting EGFR and cMET with antitumor activity.

Bafisontamab (EMB-01) induces significantly stronger anti-tumor activity in PDX and CDX tumor models than that of anti-EGFR mAb or anti-cMet mAb alone. Bafisontamab shows significant anti-tumor activity in anti-EGFR mAb resistant tumor model[1].

References:

[1]. Ren F, et al. EMB-01: An innovative bispecific antibody targeting EGFR and cMet on tumor cells mediates a novel mechanism to improve anti-tumor efficacy. Cancer Research, 2020,

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

80(16_Supplement): 528-528.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA