

Product Data Sheet

Product Name: HBF-0259

Cat. No.: GC63003

Chemical Properties

Cas. No. 957011-15-3


Formula $C_{16}H_{12}Cl_2FN_5$ M.Wt 364.2

Solubility DMSO : 250 mg/mL (686.44 mM; Storage Store at -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

HBF-0259 is a potent and selective inhibitor of hepatitis B virus (HBV) surface antigen (HBsAg) secretion, with an EC₅₀ of 1.5 μM in HepG2.2.15 cells. HBF-0259 has no effect on HBV DNA

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

synthesis[1][2].

HBF-0259 (0.0158-50.0 μ M) inhibits the HBsAg secretion (EC_{50} =1.5 μ M) in HepDE19 cells, with an CC_{50} of >50 μ M[1]. HBF-0259 (2-8 μ M) significantly reduces the abundance of both L and M in the supernatants of HepG2.2.15 cells[1].

[1]. Dougherty AM, et, al. A substituted tetrahydro-tetrazolo-pyrimidine is a specific and novel inhibitor of hepatitis B virus surface antigen secretion. Antimicrob Agents Chemother. 2007 Dec;51(12):4427-37.

[2]. Yu W, et, al. Design, synthesis, and biological evaluation of triazolo-pyrimidine derivatives as novel inhibitors of hepatitis B virus surface antigen (HBsAg) secretion. J Med Chem. 2011 Aug 25;54(16):5660-70.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA