

Product Data Sheet

Product Name: 4,6-Bis(3,5-di(pyridin-3-yl)phenyl)-2-methylpyrimidine

Cat. No.: GD05352

Chemical Properties

Cas. No. 925425-96-3

Formula C₃₇H₂₆N₆

M.Wt 554.643


Solubility

Storage 2-8°C argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA