

Product Data Sheet

Product Name: p53 Activator 7

Cat. No.: GC69648

Chemical Properties

Cas. No. 2849340-59-4


Formula $C_{27}H_{32}F_3N_4OP$ M.Wt 516.54

Solubility DMSO : 100 mg/mL (193.60 mM; Storage 4°C, away from moisture and light
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA