

Product Data Sheet

Product Name: Boc-Val-Pro-Arg-AMC

Cat. No.: GA21166

Chemical Properties

Cas. No. 65147-04-8

Formula $C_{31}H_{45}N_7O_7$

M.Wt 627.74

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Boc-VPR-AMC, a specific, highly fluorogenic substrate for thrombin and the thrombin-staphylocoagulase complex. The thrombin assay has been used for the indirect determination of coagulation factor III. Trypsin, trypsin from rat mast cells, and two try

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA