

Product Data Sheet

Product Name: CamA-IN-1

Cat. No.: GC68823

Chemical Properties

Cas. No.


Formula $C_{23}H_{36}N_6O_6$ M.Wt 492.57

Solubility DMSO : 100 mg/mL (203.02 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA