

## Product Data Sheet

Product Name: H-Phe-pNA

Cat. No.: GA10096

### Chemical Properties

Cas. No. 2360-97-6

Chemical Name (2S)-2-amino-N-(4-nitrophenyl)-3-phenylpropanamide

SMILES C1=CC=C(C=C1)CC(C(=O)NC2=CC=C(C=C2)[N+](=O)[O-])N


Formula  $C_{15}H_{15}N_3O_3$  M.Wt 285.3

Solubility  $\geq 28.5\text{mg/mL}$  in DMSO,  $\geq 19.6$  mg/mL in Water with ultrasonic Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



### Background

H-Phe-pNA

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA