

Product Data Sheet

Product Name: cis-5-Norbornene-exo-2,3-dicarboxylic anhydride

Cat. No.: GD06082

Chemical Properties

Cas. No. 2746-19-2

Formula C₉H₈O₃

M.Wt 164.158


Solubility

Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Himic anhydride is an anhydride modified cantharidin analogue that inhibits protein phosphatase 2A. Himic anhydride can be used in asthma research.

[1]. McCuskey A, et al. Anhydride modified cantharidin analogues. Is ring opening important in the

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

inhibition of protein phosphatase 2A? Eur J Med Chem. 2000 Oct;35(10):957-64.

[2]. Rosenman KD, et al. Occupational asthma caused by himic anhydride. Scand J Work Environ Health. 1987 Apr;13(2):150-4.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA