

Product Data Sheet

Product Name: PF-3644022

Cat. No.: GC14645

Chemical Properties

Cas. No. 1276121-88-0

Chemical Name (R)-10-methyl-3-(6-methylpyridin-3-yl)-11,12-dihydro-10H-[1,4]diazepino[5',6':4,5]thieno[3,2-f]quinolin-8-ol

SMILES CC1=NC=C(C2=NC(C=CC(S3)=C4C5=C3C(O)=N[C@](CN5)([H])C)=C4C=C2)C=C1

Formula $C_{21}H_{18}N_4OS$

M.Wt 374.46


Solubility <14.98mg/ml in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

PF-3644022 is a potent and selective mitogen-activated protein kinase (MAPK)-activated protein kinase-2 (MK2) inhibitor ($K_i = 3$ nM). [1]

PF-3644022 is the first oral MK2 inhibitor in both acute and chronic models of inflammation. PF-3644022 reversibly ATP-competitively inhibits MK2 enzyme activity with good selectivity across 200 human kinases. [1]

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

In the human U937 monocytic cell line or peripheral blood mononuclear cells, PF-3644022 potently inhibits TNF production with IC₅₀ value of 160 nM. In LPS-stimulated human whole blood, PF-3644022 blocks TNF and IL-6 production with IC₅₀ values of 1.6 and 10.3 μ M, respectively. In U937 cells and blood, Inhibition of TNF correlates closely with inhibition of phospho-heat shock protein 27, a target biomarker of MK2 activity. [1]

In the rat LPS-induced TNF model, PF-3644022 efficacy in the chronic inflammation model is strongly correlated with maintaining a Cmin higher than the EC₅₀ measured. PF-3644022 exhibits efficaciously oral activity and good pharmacokinetic parameters in both the chronic streptococcal cell wall-induced arthritis model and the rat acute LPS-induced TNF model. [1]

References:

1. Mourey RJ, Burnette BL, Brustkern SJ et al. A benzothiophene inhibitor of mitogen-activated protein kinase-activated protein kinase 2 inhibits tumor necrosis factor alpha production and has oral anti-inflammatory efficacy in acute and chronic models of inflammation. J Pharmacol Exp Ther. 2010 Jun;333(3):797-807. doi: 10.1124/jpet.110.166173. Epub 2010 Mar 17.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA