

Product Data Sheet

Product Name: VH 101, acid

Cat. No.: GC71570

Chemical Properties

Cas. No. 2408341-97-7

Formula $C_{28}H_{35}FN_4O_6S$


M.Wt 574.66

Solubility Storage 4°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

VH 101, acid is a functionalized von-Hippel-Lindau (VHL) protein ligand for PROTAC research and development. VH 101, acid contains an E3 ligase ligand plus an alkyl linker with terminal amine ready for conjugation to a target protein ligand.

References:

[1]. Jian Jin, et al. Protein tyrosine kinase 6 (ptk6) degradation / disruption compounds and methods of use. WO2020010204A1.2020-01-09.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA