

Product Data Sheet

Product Name: Febuxostat 67M-1

Cat. No.: GC79625

Chemical Properties

Cas. No. 887945-96-2

Formula $C_{16}H_{16}N_2O_4S$

M.Wt 332.37


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Febuxostat 67M-1 is an active metabolite of Febuxostat. Febuxostat 67M-1 retains inhibitory activity against Xanthine Oxidase, but its plasma concentration is much lower than that of Febuxostat.

Febuxostat 67M-1 can be used in the research of hyperuricemia [1].

References:

[1]. Xie H, et al. An HPLC-MS/MS method for simultaneous determination of the active metabolites of febuxostat (67M-1, 67M-2 and 67M-4) in human plasma. Journal of chromatography. B, Analytical technologies in the biomedical and life sciences. 2014 Nov 01;970:24-30.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA