

Product Data Sheet

Product Name: TH470
Cat. No.: GC73353

Chemical Properties

Cas. No. 2834739-51-2


Formula $C_{30}H_{31}N_5O_5S_2$ M.Wt 605.73

Solubility DMSO : 100 mg/mL (165.09 mM; Storage -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

TH470 is a highly selective LIMK1/2 (LIM kinase1/2) inhibitor (LIMK1 IC₅₀=9.8 nM; LIMK2 IC₅₀=13 nM), and can be used in orphan disease research.

TH470 (0.05-5 μM; 12 h) inhibits the growth of neurite[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Hanke T, et al. Development and Characterization of Type I, Type II, and Type III LIM-Kinase Chemical Probes. J Med Chem. 2022 Oct 13;65(19):13264-13287.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA