

Product Data Sheet

Product Name: (3S,4R)-GNE-6893

Cat. No.: GC73214

Chemical Properties

Cas. No. 2415374-97-7


Formula $C_{23}H_{24}FN_5O_4$ M.Wt 453.47

Solubility DMSO : 100 mg/mL (220.52 mM; Storage -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

(3S,4R)-GNE-6893 is a potent and orally active HPK1 inhibitor. (3S,4R)-GNE-6893 can be used in research of cancer.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Andreas S, et, al. A Convergent Synthesis of HPK1 Inhibitor GNE-6893 via Palladium-Catalyzed Functionalization of a Tetrasubstituted Isoquinoline. Org. Process Res. Dev. 2023;27(3):523-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA