

Product Data Sheet

Product Name: Magnesium iodide hydrate

Cat. No.: GD17926

Chemical Properties

Cas. No. 14332-62-8

Formula $MgI_2 \cdot xH_2O$

M.Wt 278.11


Solubility

Storage Store at room temperature argon filling, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA