

Product Data Sheet

Product Name: Z-Arg-OH.HCl

Cat. No.: GA10045

Chemical Properties

Cas. No. 56672-63-0

Chemical Name (2S)-5-(diaminomethylideneamino)-2-(phenylmethoxycarbonylamino)pentanoic acid;hydrochloride

SMILES C1=CC=C(C=C1)COC(=O)NC(CCCN=C(N)N)C(=O)O.Cl


Formula $C_{14}H_{20}N_4O_4 \cdot HCl$ M.Wt 344.8

Solubility Soluble in water or 1% acetic acid Storage Desiccate at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Z-Arg-OH.HCl

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA