

Product Data Sheet

Product Name: Guanosine 5'-(tetrahydrogen triphosphate) Sodium Salt

Cat. No.: GF11080

Chemical Properties

Cas. No. 24905-71-3

Formula C₁₀H₁₆N₆O₁₄P₃xNa

M.Wt 523.18

Solubility

Storage 20°C,

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA