

Product Data Sheet

Product Name: Fmoc-Asp(OtBu)-Wang resin

Cat. No.: GA10064

Chemical Properties

Cas. No. N/A

Chemical Name (S,E)-tert-butyl 3-((((9H-fluoren-9-yl)methoxy)(hydroxy)methylene)amino)-4-oxobutanoate

SMILES CC(C)(OC(C[C@@]([H])(C=O)O)N=C(OCC1C2=CC=CC=C2C3=CC=CC=C31)\O)([H])C=O)=O)C

Formula C₂₃H₂₅NO₅

M.Wt 395.45


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-Asp(OtBu)-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA