

Product Data Sheet

Product Name: H-Arg-βNA · HCl

Cat. No.: GA22099

Chemical Properties

Cas. No. 18905-73-2

Formula $C_{16}H_{21}N_5O \cdot HCl$

M.Wt 335.84

Solubility 50 mg/mL in methanol

Storage Store at -20°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Substrate for aminopeptidase B (arginine aminopeptidase) and cathepsin H.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA