

Product Data Sheet

Product Name: Thalidomide-5-O-C6-NH₂ hydrochloride

Cat. No.: GC73248

Chemical Properties

Cas. No. 2761398-05-2

Formula C₁₉H₂₄ClN₃O₅


M.Wt 409.86

Solubility Storage -20°C, protect from light, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Thalidomide-5-O-C₆-NH₂ hydrochloride is the Thalidomide-based cereblon ligand used in the recruitment of CRBN protein. Thalidomide-5-O-C₆-NH₂ hydrochloride can be connected to the ligand for protein by a linker to form PROTACs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA