

Product Data Sheet

Product Name: Propargyl ether

Cat. No.: GD20328

Chemical Properties

Cas. No. 6921-27-3

Formula C_6H_6O

M.Wt 94.11


Solubility

Storage $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Propargyl Ether is a useful ether for proteomics research. Propargyl Ether can be used for the

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

preparation of polyferrocenes^[1].

References:

[1].Crich D, et al. Enhanced diastereoselectivity in beta-mannopyranosylation through the use of sterically minimal propargyl ether protecting groups. J Org Chem. 2006 Apr 14;71(8):3064-70.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA