

Product Data Sheet

Product Name: Thalidomide-5-PEG6-NH₂ hydrochloride

Cat. No.: GC73556

Chemical Properties

Cas. No.

Formula C₂₅H₃₆ClN₃O₁₀

M.Wt 574.02


Solubility

Storage -20°C, sealed storage, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Thalidomide-5-PEG6-NH₂ drochloride is the Thalidomide -based cereblon ligand used in the recruitment of CRBN protein. Thalidomide-5-PEG6-NH₂ drochloride can be connected to the ligand for protein by a linker to form PROTACs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA