

Product Name: PROTAC GPX4 degrader-1
Cat. No.: GC73223

Cas. No. 2916433-81-1

Formula $C_{50}H_{57}ClN_{10}O_{10}$


M.Wt 993.5

Solubility	DMSO : 100 mg/mL (100.65 mM; Need ultrasonic)	Storage	4°C, stored under nitrogen
------------	---	---------	----------------------------

General tips	For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.
--------------	---

Shipping Condition	Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.
--------------------	--

Structure



PROTAC GPX4 degrader-1 (DC-2) is a PROTAC-based GPX4 degrader, with a DC50 of 0.03 μ M in HT1080 cells.

PROTAC GPX4 degrader-1 (DC-2) has been found to induce GPX4 degradation and subsequent ferroptosis[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Han Wang, et al. Discovery of ML210-Based glutathione peroxidase 4 (GPX4) degrader inducing ferroptosis of human cancer cells. Eur J Med Chem. 2023 Jun 5;254:115343.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA