

Product Data Sheet

Product Name: GW-870086

Cat. No.: GC31936

Chemical Properties

Cas. No. 827319-43-7

SMILES O=C(C1C(C)(C)C1(C)C)O[C@]2([C@@]3([C@@]([C@@]4([H])[C@@](F)([C@@]5(C([C@@H](F)C4)=CC(C=C5)=O)C)[C@@H](O)C3)([H])C[C@H]2C)C)C(OCC#N)=O

Formula $C_{31}H_{39}F_2NO_6$ M.Wt 559.64

Solubility DMSO : 100 mg/mL (178.69 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.


Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Protocol

Cell experiment: Cells are exposed to DEX, PRED, or GW-870086 at 0, 1, 3, 10, 30, 100 or 300 nM for 5 days. Compounds (GW-870086, etc.) are reconstituted in DMSO as 300-mM stock solutions that are diluted to give final drug concentrations in DMEM. Media (DMEM supplemented with 1% FBS) is changed after 48 hours and equal volumes of media from each group are collected after the final 72 hours of treatment, concentrated ($\times 10$), and added to an equal volume of $\times 2$ sample buffer. Cells are then rinsed with PBS and scraped directly into 100- μ L $\times 2$ sample buffer/well. All samples are incubated at 100°C for 8 minutes and stored at -20°C until analysis[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

- [1]. Uings IJ, et al. Discovery of GW870086: a potent anti-inflammatory steroid with a unique pharmacological profile. Br J Pharmacol. 2013 Jul;169(6):1389-403.
- [2]. Stamer WD, et al. Unique response profile of trabecular meshwork cells to the novel selective glucocorticoid receptor agonist, GW870086X. Invest Ophthalmol Vis Sci. 2013 Mar 1;54(3):2100-7.

Background

GW-870086 is a potent anti-inflammatory agent, acting as a glucocorticoid receptor agonist, with a pIC₅₀ of 10.1 in A549 cells expressing NF-κB.

GW-870086 a glucocorticoid receptor agonist, inhibits NFκB reporter gene with a pIC₅₀ of 10.1, but shows no effect on the MMTV reporter genes in A549 cells, and has little or no activity at the oestrogen receptor, progesterone receptor, progesterone receptor, mineralocorticoid receptor or androgen receptor. GW-870086 (1 pM-1 μM) dose-dependently inhibits the L-6 release induced by TNF-α in A549 epithelial carcinoma cells and by IL-1 in MG63 osteosarcoma cells (pIC₅₀s, 9.6, 10.2) [1]. GW-870086 (GW870086X; 10-100 nM) significantly increases fibronectin secretion. However, GW-870086 has no effect on MMP2 secretion, and does not increase cellular myocilin[2].

[1]. Uings IJ, et al. Discovery of GW870086: a potent anti-inflammatory steroid with a unique pharmacological profile. Br J Pharmacol. 2013 Jul;169(6):1389-403. [2]. Stamer WD, et al. Unique response profile of trabecular meshwork cells to the novel selective glucocorticoid receptor agonist, GW870086X. Invest Ophthalmol Vis Sci. 2013 Mar 1;54(3):2100-7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA