

Product Data Sheet

Product Name: PHI-101
Cat. No.: GC73497

Chemical Properties

Cas. No. 2127107-15-5


Formula $C_{19}H_{19}FN_4O_2S$ M.Wt 386.44

Solubility DMSO : 100 mg/mL (258.77 mM; Need ultrasonic) Storage $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PHI-101 is an orally active FLT3 inhibitor that overcomes resistance to multiple drug-resistant mutations. PHI-101 potently inhibits FLT3 single activating mutations (ITD or TKD mutants) and has inhibitory activity against FLT3 double (ITD/D835Y or ITD/F691L) and triple (ITD/D835Y/F691L)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

resistance mutations. PHI-101 has potential for research in relapsed or refractory acute myeloid leukemia (AML).

References:

[1]. Nam K Y, et al. PHI-101 Is a Potent Third-Generation FLT3 Inhibitor Developed to Overcome Resistance in Acute Myeloid Leukemia[J]. Blood, 2020, 136: 28.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA