

Product Data Sheet

Product Name: AKI603
Cat. No.: GC62433

Chemical Properties

Cas. No. 1432515-73-5


Formula $C_{19}H_{23}N_9O_2$ M.Wt 409.45

Solubility DMSO : 125 mg/mL (305.29 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

AKI603 is an inhibitor of Aurora kinase A (AurA), with an IC₅₀ of 12.3 nM. AKI603 is developed to

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

overcome resistance mediated by BCR-ABL-T315I mutation. AKI603 exhibits strong anti-proliferative activity in leukemic cells[1][2].

AKI603 (0.039-0.6 μ M; 48 hours) extensively inhibits proliferation of leukemia cells[1].AKI603 (0.039-0.6 μ M; 48 hours) significantly inhibits the phosphorylation of AurA in NB4, K562, and Jurkat cell lines in a dose-dependent manner while the level of total AurA protein is not changed[1].AKI603 inhibits the proliferation and colony formation of imatinib resistant CML cells[1].AKI603 (0.3-0.6 μ M; 48 hours) inhibits cell proliferation and colony formation capacities in imatinib-resistant CML cells by inducing cell cycle arrest with polyploidy accumulation[1].Inhibition of AurA by AKI603 induces leukemia cell senescence in both BCR-ABL wild type and T315I mutation cells[1].AKI603 exhibits inhibitory activities on breast cancer cell proliferation, such as SUM149 (IC₅₀=2.04), BT549 (IC₅₀=0.86), MCF-7 (IC₅₀=0.97), MCF-7-Epi (IC₅₀=21.01), Sk-br-3 (IC₅₀=0.73), MDA-MB-231 (IC₅₀=3.49), MDA-MB-453 (MTT, IC₅₀=0.18; Cell counting, IC₅₀=0.19), MDA-MB-468 (MTT, IC₅₀=0.15; Cell counting, IC₅₀=0.17)[2].

AKI603 (12.5-25 mg/kg; i.p.; every 2 days; for 14 days) abrogates the growth of xenografted KBM5-T315I cells in nude mice[1].AKI603 exhibits moderate oral bioavailability (rat 28.7%) and C_{max} (rat 202.4 μ g/L) following oral administration (rat 25 mg/kg)[3].AKI603 exhibits terminal elimination half-life (rat 8.9 h) following intravenous administration (rat 2.5 mg/kg)[3].

[1]. Le-Xun Wang, et al. Aurora A Kinase Inhibitor AKI603 Induces Cellular Senescence in Chronic Myeloid Leukemia Cells Harboring T315I Mutation. Sci Rep. 2016 Nov 8;6:35533.

[2]. Fei-Meng Zheng, et al. A novel small molecule aurora kinase inhibitor attenuates breast tumor-initiating cells and overcomes drug resistance. Mol Cancer Ther. 2014 Aug;13(8):1991-2003.

[3]. Zhenzhen Zhao, et al. Determination of a novel Aurora-A (AurA) kinase AKI603 by UPLC-MS/MS and its application to a bioavailability study in rat. J Pharm Biomed Anal. 2016 Jun 5;125:303-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA