

Product Data Sheet

Product Name: (+)-Medioresinol

Cat. No.: GC67931

Chemical Properties

Cas. No. 40957-99-1

Formula $C_{21}H_{24}O_7$

M.Wt 388.41


Solubility DMSO : 50 mg/mL (128.73 mM;
ultrasonic and warming and heat to
60°C)

Storage 4°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA