

Product Data Sheet

Product Name: PI3K/mTOR ligand-1

Cat. No.: GC81042

Chemical Properties

Cas. No. 1197160-66-9


Formula $C_{25}H_{27}N_7O_5$ M.Wt 505.53

Solubility DMSO: 25 mg/mL (49.45 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

PI3K/mTOR ligand-1 is a ligand for the target protein of PROTAC (PI3K / mTOR). PI3K/mTOR ligand-1 can be used to synthesize GP262 [1].

In Vitro, PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

References:

[1]. Zhu W, et al. Discovery of the First Potent PI3K/mTOR Dual-Targeting PROTAC Degradator for Efficient Modulation of the PI3K/AKT/mTOR Pathway. J Med Chem. 2025 Dec 25;68(24):26188-26205.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA