

Product Data Sheet

Product Name: (2S, 4R)-4-Phenylpyrrolidine-2-carboxylic acid hydrochloride

Cat. No.: GF52883

Chemical Properties

Cas. No. 82087-67-0

Formula C₁₁H₁₄ClNO₂

M.Wt 227.69

Solubility

Storage 2-8°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA