

Product Data Sheet

Product Name: Fmoc-4-methoxy-L-homophenylalanine

Cat. No.: GF11440

Chemical Properties

Cas. No. 205182-52-1

Formula C₂₆H₂₅NO₅

M.Wt 431.4804

Solubility

Storage Store at RT

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA