

Product Data Sheet

Product Name: DOTAP mesylate

Cat. No.: GB58331

Chemical Properties

Cas. No. 252769-92-9


Formula C₄₃H₈₃NO₇S M.Wt 758.19

Solubility Ethanol : 33.33 mg/mL (43.96 mM; Storage
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA