

Product Data Sheet

Product Name: CTR/AMYR modulator-2

Cat. No.: GC79274

Chemical Properties

Cas. No. 3108095-52-6

Formula $C_{33}H_{35}ClF_2N_8O_3$

M.Wt 665.13

Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

CTR/AMYR modulator-2 is a dual agonist of human amylin receptor 3 (AMY3R) and calcitonin receptor (CTR). CTR/AMYR modulator-2 activates AMY3R and CTR, stimulating Rel EC 50 of cAMP production in UMUC3 cells by 4 nM and 1 nM, respectively. CTR/AMYR modulator-2 may be used in research on type 2 diabetes and obesity ^[1].

In Vitro, CTR/AMYR modulator-2 (Example 107) acts as an agonist at the human amylin receptor (AMY3R) in UMUC3 cells, with a Rel EC₅₀ of 4 nM and an E_{max} of 99% as determined in a 20-point concentration-response assay following 30 min of incubation^[1]. CTR/AMYR modulator-2 acts as an

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

agonist at the human calcitonin receptor (CTR) in UMUC3 cells, with a Rel EC₅₀ of 1 nM and an E_{max} of 88% under the same 20-point assay conditions^[1].

References:

[1]. Renato A. BAUER, et al. 4-phenyl-1h-pyrazolo[3,4-b]pyridine-5-carboxamide derivatives as amylin and/or calcitonin receptor agonists for the treatment of type 2 diabetes, obesity or overweight. WO2025240893A1. 2025.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA