

Product Data Sheet

Product Name: PI4K-IN-1

Cat. No.: GC72785

Chemical Properties

Cas. No. 1800017-49-5


Formula $C_{24}H_{27}N_3O_3S$ M.Wt 437.55

Solubility DMSO : 125 mg/mL (285.68 mM; Storage -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PI4K-IN-1 (compound 44) is a potent PI4KIII inhibitor, with pIC50 values of 9.0 and 6.6 for PI4KIII α and PI4KIII β , respectively. PI4K-IN-1 also inhibits PI3K $\alpha/\beta/\gamma/\delta$, with pIC50 values of

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

4.0/<3.7/5.0/<4.1, respectively.

References:

[1]. Raubo P, et al. Discovery of potent, selective small molecule inhibitors of α -subtype of type III phosphatidylinositol-4-kinase (PI4KIII α). Bioorg Med Chem Lett. 2015 Aug 15;25(16):3189-93.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA