

Product Data Sheet

Product Name: (R)-BPO-27

Cat. No.: GC31665

Chemical Properties

Cas. No. 1415390-47-4

SMILES O=C(C1=CC=C(O[C@@H](C2=CC=C(Br)O2)C3=C(N(C)C4=O)C(C(N4C)=O)=C(C5=C C=CC=C5)N36)C6=C1)O


Formula $C_{26}H_{18}BrN_3O_6$ M.Wt 548.34

Solubility DMSO : ≥ 14.28 mg/mL (26.04 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Protocol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Cell experiment: Whole-cell recordings are done on CFTR-expressing CHO-K1 cells. After establishing the whole-cell configuration, BPO-27 is added for 5 minutes, and then CFTR is activated by the addition of forskolin (10 μ M) in the continued presence of BPO-27 (0.5 or 1 μ M). Whole-cell currents are elicited by applying hyperpolarizing and depolarizing voltage pulses from a holding potential of 0 mV to potentials between +80 and -80 mV in steps of 20 mV. Recordings are made at room temperature using an Axopatch-200B. Currents are digitized with a Digidata 1440A converter and filtered at 5 kHz[2].

Animal experiment: Rats: (R)-BPO-27 is formulated at 1 mg/mL in 5% DMSO, 2.5% Tween-80 and 2.5% PEG400 in water. Male mice in a CD1 genetic background are administered 300 μ L of the (R)-BPO-27 formulation by intraperitoneal injection. At specified times, blood samples are collected by eye bleed. At 4 h, kidneys are removed following renal arterial perfusion with PBS. Kidneys are weighed, mixed with acetic acid and homogenized for analysis[1].

References:

- [1]. Snyder DS, et al. Absolute Configuration And Biological Properties of Enantiomers of CFTR Inhibitor BPO-27. ACS Med Chem Lett. 2013 May 9;4(5):456-459.
[2]. Kim Y, et al. Benzopyrimido-pyrrolo-oxazine-dione (R)-BPO-27 Inhibits CFTR Chloride Channel Gating by Competition with ATP. Mol Pharmacol. 2015 Oct;88(4):689-96.

Background

(R)-BPO-27 is a potent CFTR inhibitor with an IC₅₀ of 4 nM.

The benzopyrimido-pyrrolo-oxazinedione BPO-27 is an analogue of PPQ-102, which inhibits CFTR with an IC₅₀ of 8 nM. The R enantiomer of BPO-27 inhibits CFTR chloride conductance with an IC₅₀ of 4 nM, while S enantiomer is inactive. In vitro metabolic stability in hepatic microsomes shows both enantiomers as stable, with less than 5% metabolism in 4 h[1]. (R)-BPO-27 binds near the canonical ATP binding site. Whole-cell patch-clamp studies shows linear CFTR currents with a voltage-

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

independent (R)-BPO-27 block mechanism. At a concentration of (R)-BPO-27 that inhibits CFTR chloride current by 50%, the EC₅₀ for ATP activation of CFTR increases from 0.27 to 1.77 mM[2].

Following bolus interperitoneal administration in mice, serum (R)-1 decays with $t_{1/2} \approx 1.6$ h and gives sustained therapeutic concentrations in kidney[1].

[1]. Snyder DS, et al. Absolute Configuration And Biological Properties of Enantiomers of CFTR Inhibitor BPO-27. ACS Med Chem Lett. 2013 May 9;4(5):456-459. [2]. Kim Y, et al. Benzopyrimido-pyrrolo-oxazine-dione (R)-BPO-27 Inhibits CFTR Chloride Channel Gating by Competition with ATP. Mol Pharmacol. 2015 Oct;88(4):689-96.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA