

Product Data Sheet

Product Name: BI-1408
Cat. No.: GC70338

Chemical Properties

Cas. No. 2231075-94-6


Formula $C_{22}H_{23}FN_6$ M.Wt 390.46

Solubility Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available sizes: ship with RT, or blue ice upon request.

Structure



Background

BI-1408 is a potent γ secretase modulator with an IC_{50} of 0.04 μ M for A β 42.

References:

[1]. Gerlach K, et al. Discovery of tetrahydroindazoles as a novel class of potent and in vivo efficacious gamma secretase modulators. Bioorg Med Chem. 2018 Jul 23;26(12):3227-3241.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA