

Product Data Sheet

Product Name: 5-F-UTP trisodium

Cat. No.: GC79096

Chemical Properties

Cas. No.

Formula $C_9H_{11}FN_2Na_3O_{15}P_3$ M.Wt 568.08

Solubility H_2O : 50 mg/mL (88.02 mM; Need ultrasonic) Storage Store at $-20^{\circ}C$, away from moisture

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

5-F-UTP (5-Fluorouridine-5'-triphosphate) trisodium is a monomeric raw material that can be used for nucleic acid synthesis.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA