

Product Data Sheet

Product Name: Thalidomide-4-O-C14-NH₂ hydrochloride

Cat. No.: GC73259

Chemical Properties

Cas. No.

Formula C₂₇H₄₀ClN₃O₅

M.Wt 522.08


Solubility

Storage 4°C, protect from light, stored under nitrogen

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Thalidomide-4-O-C14-NH₂ drochloride is the Thalidomide -based cereblon ligand used in the

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

recruitment of CRBN protein. Thalidomide-4-O-C14-NH₂ drochloride can be connected to the ligand for protein by a linker to form PROTACs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA