

Product Data Sheet

Product Name: WH-4-025

Cat. No.: GC62142

Chemical Properties

Cas. No. 1876463-35-2


Formula $C_{39}H_{38}F_3N_7O_5$ M.Wt 741.76

Solubility DMSO : 125 mg/mL (168.52 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

WH-4-025 is a Salt-inducible kinase (SIK) inhibitor (WO2016023014 A2)[1].

[1]. Brookline Avenue, et al. USES OF SALT-INDUCIBLE KINASE (SIK) INHIBITORS. Patent WO2016023014A2.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA