

Product Data Sheet

Product Name: Protein Kinase P34 (cd2) Substrate

Cat. No.: GA23390

Chemical Properties

Cas. No. 135546-44-0

Formula $C_{106}H_{172}N_{32}O_{32}$

M.Wt 2406.73

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

CSH-103, ADAQHATPPKKRKVEDPKDF, based on the simian virus 40 large T antigen sequence, is a specific substrate for p34 cdc2 protein kinase ($K_m = 0.07$ mM). It is not phosphorylated by other kinases (casein kinase I/II, cAMP-dependent protein kinase, glyc

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA