

Product Data Sheet

Product Name: Fmoc-D-isoGln-OH

Cat. No.: GA11054

Chemical Properties

Cas. No. 292150-20-0

Chemical Name (R,E)-4-((((9H-fluoren-9-yl)methoxy)(hydroxy)methylene)amino)-5-hydroxy-5-iminopentanoic acid

SMILES N=C(O)[C@@]([N=C(OCC1C2=CC=CC=C2C3=CC=CC=C31)\O])([H])CCC(O)=O

Formula $C_{20}H_{20}N_2O_5$

M.Wt 368.38


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-D-isoGln-OH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA