

Product Data Sheet

Product Name: Head activator

Cat. No.: GA22322

Chemical Properties

Cas. No. 79943-68-3

Formula $C_{54}H_{84}N_{12}O_{14}$

M.Wt 1125.33

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Morphogenic peptide from hydra. The neuropeptide head activator is a high-affinity ligand for the orphan G-protein-coupled receptor GPR37.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA