

Product Data Sheet

Product Name: Boc-Leu-Gly-Arg-AMC

Cat. No.: GA21044

Chemical Properties

Cas. No. 65147-09-3

Formula $C_{29}H_{43}N_7O_7$ M.Wt 601.7

Solubility DMSO : 100 mg/mL (166.20 mM; Need ultrasonic); H_2O : < 0.1 mg/mL (insoluble) Storage -20°C, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Boc-LGR-AMC, a specific fluorogenic substrate for complement component C3/C5 convertases, coagulation factor Xa, and soybean trypsin-like enzyme. It is also hydrolyzed by macropain, a high molecular weight thiol proteinase from human erythrocytes. Stock

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA