

Product Data Sheet

Product Name: m-PEG3-Amine

Cat. No.: GC39213

Chemical Properties

Cas. No. 74654-07-2

SMILES COCCOCCOCCN


Formula $C_7H_{17}NO_3$ M.Wt 163.21

Solubility Soluble in DMSO Storage 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

m-PEG3-Amine is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs) [1].

[1]. Nianhe Han, et al. The derivant of aplysiatoxin 10 and application thereof. CN106279352A.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA